

CRITICAL CARE DEVICE

TOURNIQUET | DIATHERMY | VENTILATOR





SURGICAL DIATHERMY



Features

- Monopolar cut: max 400W, Low, Pure, Blend & Endo Cut
- Monopolar coag: 120W, Desiccate, LCF Fulgurate, HCF Fulgurate & Spray
- Bipolar coag: 120W, Macro, Standard and Precise
- Bipolar Cut: 100W
- Vessel Seal : L-6, 150W, Auto-stop & regrasp Indicators
- Programmable power settings: 10 programs
- Tissue sensing technology
- Patient plate contact quality monitoring
- HF leakage monitoring

Specifications

Monopolar Cut	In 300 Ω
Blend	200 W
Pure	400 W
Low	300 W
Endo cut	200 W
Monopolar Coagulation	In 500 Ω
Spray	120 W
Fulgurate LCF	120 W
Fulgurate HCF	120 W
Desiccate	120 W
Bipolar Coagulation	In 100 Ω
Macro	120 W
Standard	070 W
Precise	120 W
Bipolar Cut	In 100 Ω
Pure Cut	120 W
Vessel Sealer	In 100 Ω
Levels	6 Levels
Power	Up to 150W



Specifications

Monopolar Cut	In 300 Ω
Blend	200 W
Pure	400 W
Low	300 W
Endo cut	200 W
Monopolar Coagulation	In 500 Ω
Spray	120 W
Fulgurate LCF	120 W
Fulgurate HCF	120 W
Desiccate	120 W
Bipolar Coagulation	In 100 Ω
Macro	120 W
Standard	070 W
Precise	120 W
Bipolar Cut	In 100 Ω
Pure Cut	120 W

Features

- Monopolar cut: max 400W, Low, Pure, Blend & Endo
- Monopolar coag: 120W, Desiccate, LCF Fulgurate, HCF Fulgurate & Spray
- Bipolar coag: 120W, Macro, Standard and Precise
- Bipolar Cut: 100W
- Programmable power settings: 10 programs
- Tissue sensing technology
- Patient plate contact quality monitoring
- HF leakage monitoring





Specifications

Monopolar Cut	In 300 Ω
Blend	200 W
Pure	300 W
Monopolar Coagulation	In 500 Ω
Spray	100 W
Fulgurate	100 W
Desiccate	120 W
Bipolar Coagulation	In 125 Ω
Standard	070 W
Bipolar Cut	In 300 Ω
Pure Cut	120 W

Features

- Monopolar cut: max 300W, 3 different blends.
- Monopolar coagulation: max 100W, Desiccate, Fulgurate & Spray modes
- Bipolar coagulation: 70W.
- Bipolar Cut: 120 W
- Programmable power settings: 10 Programs
- Patient plate contact quality monitoring

ACCESSORIES

Silicon Split Patient Plate	1	1	1
Cable for patient plate	1	1	1
Disposable handswitch, 3pin	1	1	1
Electrode set	1	1	1
Bipolar Forcep	1	1	1
Bipolar Cable	1	1	1
Monopolar foot switch	1	-	-
Bipolar Foot switch	1	-	-
Three pedal Foot switch	-	1	1
Vessel sealing laparoscopic instrument with cable	-	-	1
Biclip With Cable	-	-	1
Power Cord	1	1	1

Tailor 300

Tailor 400

Tailor Pro



FEATURES

- SPO₂ interface (All modes), Nebulizer sync (All modes with Period/continuous function), Ultrasonic O₂ sensor
- O₂ input pressure display
- Variable i-time in S/T mode based on expiration, CPAP with expiratory relief (3 levels with patient sync)
- Alarm limits for O₂, Volume, Pressure, SPO₂, Pulse, O₂ input pressure
- Open circuit detection saves O₂ consumption
- S/T mode with Auto **start** on/off function
- Optional Invasive block with volume & pressure SIMV- V, SIMV-P, PSV Modes
- Trend & log for 48 hrs
- Proximal pressure measurement for precise pressure control
- Low flow O₂ input (0 to 15 lpm) for O₂ therapy
- Interfaced directly with any Saviour ventilator
- Interfacing Saviour needs only one Pneumatic and one electrical connection
- Operational while transportation for emergency purpose such as in ambulance
- Precise electronic peep control
- Compressed gas line not required
- Flow monitoring range +/- 180 lpm
- Volume delivery 2ml - 2000 ml

SPECIFICATIONS

- **Non Invasive Modes (6)** : CPAP, O₂ therapy, S/T, S, T, PRVC (Nebulizer sync in all modes)
- **Volume (ml)**: Adult: 200 to 2000 ml (PRVC mode)
- **Rate**: 3 to 60 bpm
- **Trigger level**: 2 to 8 lpm (S/T, S mode)
- **Inspiration flow** : 200 lpm (Total leak in patient circuit)
- **Pressure** : 4 to 45 cm H₂O
- **I-Time**: 0.3 to 3 Secs
- **O₂ supply**: 0.0 to 160 lpm (All modes with blending)
- **Monitoring Parameters**: Pip, Peep, Volume (Vt), Rate, Minute volume, leak, O₂, SPO₂, Pulse.
- **Battery Backup**: 1.5 Hrs
- **Screen size**: 7"
- **Invasive modes** :
VCV (Volume control with assist), PCV (Pressure control with Assist), (Nebulizer sync) SIMV-V, SIMV-P, PSV
- **Volume Control** : Adult :150 to 2000 ml
- **Pediatric** : 20 to 300 ml
- **Neonatal** : 2 to 200 ml
- **Pressure Control** : 4 to 60 cmH₂O
- **Rate**: 3 to 60 bpm
- **Inspiration flow** : 200 lpm
- **i-Time**: 0.3 to 3 Secs
- **Trigger level**: Flow : 1 to 12 lpm
Pressure : -1 to -10 cmH₂O
- **Peep** : 0 to 20 cmH₂O
- **O₂ supply**: upto 160 L/min (21 to 100%)
- **Safety release valve**: 75 cmH₂O (+/- 5)

ELECTRONIC TOURNIQUET

Staan
MADE FOR THE WORLD



ST-ET-002
DUAL PORT

- Microprocessor controlled
- Advanced leak detection alarm
- Low battery indicators
- Set time alarm
- Online pressure manipulation
- Auto and Manual operation

- Noiseless operation
- No external high pressure required
- Acoustic alarms
- Soft touch control
- Immediate pressurization of cuffs
- Pressure retained in case of power failure

ST-ET-001 SINGLE PORT



SPECIFICATION

Class & Type	Class I type B
Cuff pressure	50mm Hg - 650mm Hg
Pressure accuracy	± 10mm Hg
Timer	2 - 200 minutes
Power supply	220 Volts 50 Hz
Initial inflate time	< 20 seconds
Noise	< 55 db
Power	25W
Battery	Lithium ion 4hrs
Number of Cuffs	3 (standard)



Cuffs

Available size
8" 12" 15" 18" 24"
30" 34" 42" 44"

Stand & Basket For Mobility And Easy
Storage (Optional)



Staan Bio-med engineering private limited

#190 A, Bharathiar Road, Ganapathy, Coimbatore - 641 006, Tamil Nadu, India.

Ph : +91-422-2533 806, +91-422-2531 008 | Mob : +91-98422 19018,+91-93631 21030

email : info@staan.in | website : www.staan.in

we reserve the right to make alterations without prior notice due to continuous R & D.Minor deviations in color mentioned are possible due to technical reasons.